


With ctDNA Tumor Fraction, the Result is Clear

More confidence in positive and negative results


Most liquid biopsy tests offer a blurred image. True positive results stand out, but negative results are harder to interpret.


The clinical trial assay based on FoundationOne® Liquid CDx (IUO) can offer a clearer picture with ctDNA tumor fraction to support trial enrollment—even for wild-type biomarkers that rely on a confident negative result.

Ready for **trial enrollment for wild-type criteria**, such as *EGFR/ALK*-negative for immunotherapy



ctDNA tumor fraction is the **key driver** of liquid biopsy concordance for **tissue-detected driver alterations**.¹



 In paired tissue samples with a **negative** liquid biopsy result, **no driver mutations were seen** when **ctDNA tumor fraction ≥1%** (0/24)—offering **confidence in a true negative**.¹

- + Driver-positive** liquid biopsy can **guide trial enrollment for therapy**.
- Driver-negative** can:
 - Informative negative with ctDNA tumor fraction ≥ 1%.**
Guide trial enrollment for therapy—even for wild-type biomarkers.
 - Indeterminate negative with low ctDNA tumor fraction.**
Suggests potential for benefit with reflex tissue CGP testing where driver alterations are frequently detected.

Insights on **molecular response** add **value** and **complement imaging**.


Standard imaging to assess treatment response can take **months**, and the results can be **inconclusive**. Is a small visual change in tumor volume a signal of response, stability or progression?

FoundationOne® Monitor can add to this signal by showing **changes in ctDNA tumor fraction** to provide insights on molecular response that **complement imaging**.²



Changes in ctDNA tumor fraction from **pre-treatment** to an **early on-treatment timepoint** can provide insights into molecular response.³

Example: Retrospective analysis of mCRPC subjects showed clear stratification when ctDNA tumor fraction decreased by $\geq 75\%$ or was not detected from pre-treatment to on-treatment



ctDNA tumor fraction provides a clearer picture compared to variant allele frequency (VAF) by combining a multi-omic assessment of DNA with filtering for clonal hematopoietic (CH) variants, which can confound VAF alone.^{4,5}



CH variants can confound ctDNA assessment based on VAF alone, with examples of variants with an increase, decrease or no change in VAF across treatment timepoints—**leading to less clear picture of molecular response**.



ctDNA tumor fraction incorporates CH filtering to provide a **more precise estimate of ctDNA than VAF assessment alone**.⁶

FAST ANSWERS

Get insights on ctDNA tumor fraction quickly to add **value** to your **clinical programs** or **translational research**.

8 days
PROSPECTIVE RESULTS

14 days
RETROSPECTIVE RESULTS

The studies above used the clinical trial assay based on the Foundation Medicine liquid platform. The ctDNA tumor fraction algorithm used in these studies is for research use only.

References

1. Rolfo CD et al. Utility of ctDNA Tumor Fraction to Inform Negative Liquid Biopsy (LBx) Results and Need for Tissue Reflex in Advanced Non-Small Cell Lung Cancer (aNSCLC). J Clin Onc. 2023; 41(16): 9076-9076. DOI: 10.1200/JCO.2023.41.16_suppl.9076.
2. For research use only. Not for use in diagnostic procedures.
3. Sweeney C, Xu C, He J, et al. Evaluation of Circulating Tumor DNA Kinetics as a Prognostic Biomarker for Overall Survival in Metastatic Castrate Resistant Prostate Cancer. Cancer Res 2023; 83(7_Supplement): 3362. DOI:https://doi.org/10.1158/1538-7445.AM2023-3362.
4. Chan HT, Nagayama S, Chin YM, et al. Clinical Significance of Clonal Hematopoiesis in the Interpretation of blood Liquid Biopsy. Mol Onc. 2020; 14(8): 1719-1730. https://doi.org/10.1002/1878-0261.12727.
5. Husain H, Pavlick DC, Fendler BJ, et al. Tumor Fraction Correlates with Detection of Actionable Variants Across >23,000 Circulating Tumor DNA Samples. JCO Precis Oncol. 2022 Oct;6:e2200261. doi: 10.1200/PO.22.00261.
6. Rolfo CD, Madison R, Pasquina LW. Utility of ctDNA Tumor Fraction to Inform Negative Liquid Biopsy (LBx) Results and Need for Tissue Reflex in Advanced Non-Small Cell Lung Cancer (aNSCLC). J Clin Onc. 2023; 41(16): 9076-9076. DOI: 10.1200/JCO.2023.41.16_suppl.9076.